

BISHOP LAB DMC1 CHIP PROTOCOL

Note: This protocol is optimized for meiotic yeast cells sporulating in liquid medium.

CROSS-LINKING

1. Harvest 100mLs of cells from sporulating culture (only half of these cells will be cross-linked; the other half is reserved for the no cross-linking control – to process these cells, skip cross-linking steps and start at step 7)
2. To cross-link, aliquot 50mLs of cells into an Erlenmeyer flask
3. Add formaldehyde to 1%
4. Shake gently at room temperature for 15 minutes (swirl flask occasionally because the yeast will settle to the bottom)
5. To quench cross-linking, add glycine to final concentration of 125mM
6. Shake gently at room temperature for 5 minutes
7. Pour cross-linked or non-cross-linked cells into 50mL conical tubes
8. Spin down 5 min at 3000 rpm in tabletop centrifuge
9. Pour off supernatant and resuspend cells in 1x PBS
10. Repeat spin
11. Pour off supernatant and resuspend cell pellet in 1-2 mLs 1x PBS
12. Divide each sample evenly into 2 microcentrifuge tubes
13. Quick spin, aspirate supernatant
14. Freeze pellet on dry ice
15. Store at -70°C

Yield: 2 cross-linked pellets and 2 non-cross-linked pellets per strain

LYSIS AND SONICATION

1. Thaw cell pellet on ice for ~15 minutes
2. Add protease inhibitors (PI) fresh to 1M NaCl Lysis Buffer
3. Resuspend thawed cell pellets in 500 μL buffer + PI
4. Add cold glass beads to 500 μL mark on tube
5. Vortex at top speed for 5-10 minutes at 4°C
6. Separate cell lysate from glass beads - poke a hole in the bottom of the tube using a hot 20-gauge needle, place tube over a fresh microcentrifuge tube, and spin in for 1 minute at 3000 rpm at 4°C
7. Sonicate lysate to final shear size of 1Kb (on a Branson Sonifier 450: output control=1, duty cycle=50%, count 20 pulses. Repeat this 3 times for each sample, with at least 1 minute on ice in between sonications)
8. Spin sonicated lysate for 10 minutes at top speed in a microcentrifuge at 4°C
9. Transfer supernatant (lysate) to a fresh tube, noting volume – store on ice.
10. If necessary to increase amount of protein, resuspend cell pellet from step 8 in 500 μL 1M NaCl lysis buffer + PI and repeat steps 3-9 (In step 9, combine lysate with lysate from first round, noting final volume)
11. Measure concentration of protein in lysate using Bradford Assay (average yield is about 5 $\mu\text{g}/\mu\text{L}$ protein per sample)
12. Reserve 30-40 μL lysate for input DNA sample (2 rounds of phenol/chloroform extractions followed by EtOH precipitation of DNA)

IMMUNOPRECIPITATION

(Note: These conditions are optimized for 2mg input protein for IP with Guinea Pig Dmc1 Antibody)

1. For each sample, calculate amount of lysate for 2mg protein
2. Aliquot this volume of lysate into microcentrifuge tube
3. If volume is less than 500 μ L, add 1M NaCl lysis buffer to adjust volume to 500 μ L
4. Add 1 μ L Dmc1 Guinea Pig antibody
5. Rotate samples overnight at 4°C
6. Add 25 μ L Protein G magnetic beads (washed 2 times 1M NaCl lysis buffer) and 0.6 μ L BSA (10mg/mL stock)
7. Rotate samples for 2 hours at 4°C

WASHING BEADS AND COLLECTING CHIP EXTRACT

1. Separate beads from lysate on magnet (this will take 1-2 minutes)
2. Remove supernatant from beads, discard
3. Resuspend beads in 1mL 1M NaCl lysis buffer
4. Rotate tubes at room temperature for 5 minutes (This is critical! Skimping on wash time yields dirty IPs!)
5. Repeat wash 4 more times with 1M NaCl lysis buffer
6. Repeat wash 3 times with 0.5M NaCl Buffer
7. Repeat wash 1 time with LiCl Wash Solution
8. Repeat wash 1 time with 1x TE, pH 8.0
9. After final wash, remove all traces of wash solution from the beads
10. Resuspend beads in 100 μ L 1x TE, pH 8.0
11. To reverse cross-links, heat samples at 65°C overnight (at minimum, samples must be at 65°C for six hours to reverse cross-links)
12. Transfer each sample carefully to a PCR tube
13. Heat to 95°C for 30 minutes in PCR machine (keeps better constant temperature than a water bath)
14. Touch spin tubes
15. Separate supernatant from beads on magnet for at least 2 minutes
16. Transfer supernatant (ChIP extract) to fresh tube, being careful to not transfer beads
17. Sample is now ready for analysis
18. For a standard PCR reaction, use 2 μ L ChIP extract

SOLUTIONS

1M NaCl Lysis Buffer

50mM HEPES/KOH, pH 7.5

1M NaCl

5mM EDTA

1% Triton X-100

0.1% Na Deoxycholate

For 10mLs 1M NaCl lysis buffer, add following protease inhibitors fresh before performing lysis step:

10 μ L 1mg/mL Leupeptin (in H₂O)

50 μ L 0.2M PMSF (in EtOH)

10 μ L 1mg/mL Pepstatin (in MeOH)

Alternatively, add 100x Protease Inhibitor cocktail to final concentration of 1x

0.5M NaCl Buffer

50mM HEPES/KOH, pH 7.5

0.5M NaCl

5mM EDTA

1% Triton X-100

0.1% Na Deoxycholate

LiCl Wash Buffer

10mM Tris-HCl, pH 8

250mM LiCl

0.5% NP-40

0.5% Na Deoxycholate

5mM EDTA