

BISHOP LAB (Modified from Red Book Protocol)

YEAST TRANSFORMATION

This protocol is optimized for 8 transformations. Adjust accordingly.

PREPARING COMPETENT YEAST

1. On day before transformations, set up overnight culture by inoculating a single colony from a fresh plate into 5mL YPDA
2. A few hours before doing transformations, inoculate 195mL YPDA with entire 5mL overnight culture (i.e. 1/40 dilution)
3. Grow cells at 30°C until OD₆₀₀ 0.1-0.2 – this will take 2-4 hours depending on how your specific strain is growing
Note: Overgrowth will significantly decrease transformation efficiency!
4. Spin down cells at 3000rpm for 3 min at room temperature
5. Resuspend pellet in 50mL sterile H₂O and transfer cells to sterile 50mL conical tube
6. Spin as in step 4
7. Resuspend pellet in 5mL 0.1M LiAc
8. Spin as in step 4
9. Resuspend pellet in 400μL 0.1M LiAc (enough for 8 transformations)

TRANSFORMATIONS

1. Aliquot 50μL yeast per transformation tube
2. Add 50μL 2mg/mL carrier DNA (salmon sperm or similar)
3. Add transformation DNA and mix well by pipeting up and down
Note: Set up a negative control tube using only water.
4. Add 600μL PEG/LiAc and mix well by pipeting up and down
Note: PEG mix is viscous so be sure to pipet slowly in order to get all of PEG into the tube and well mixed with the yeast.
5. Incubate at 30°C for 30 minutes (water bath is preferable)
6. Add 70μL DMSO
7. Heat shock immediately by incubating at 42°C for 15 minutes (water bath is preferable)
Note: Set a timer! If incubated for much longer than 15 minutes at 42°C, yeast will die from heat shock.
8. Spin tubes at maximum speed in a microcentrifuge for 1 minute at room temperature
9. Resuspend in 100-400μL sterile H₂O
10. Plate 100-200μL cells per selective plate
11. Incubate at 30°C for two days

YPDA

10g yeast extract
20g bacto peptone
10mL 0.1% adenine

Dissolve in 1L H₂O
Adjust pH to 7 with NaOH
Autoclave
Add 50mL sterile 40% glucose per liter

PEG/LiAc Mix

77.4mL 50% PEG (3350)
11mL 1M LiAc
11.6mL ddH₂O